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
1. BSI. (1969) Specification for respirators for protection against harmful dusts and gases. BS2091:1969. London: British Standard Institution.
2. AS/NZS. (1994) Respiratory protective devices. AS/NZS 1715. New South Wales, Australia: Standards Australia.
3. CFR (Code of Federal Regulations). (2002) Approval of respiratory protective devices. Title 42, Part 84. Washington, DC: US Government Printing Office, Office of the Federal Register.
4. HSE. (2003) HSE operational circular OC 282/28 fit testing of respiratory protective equipment.
5. da Roza RA, Biermann AH, Foote KL *et al.* (1991) Evaluation of Portacount for determining respirator fit factors, Part III: human subject tests and comparison with an aerosol photometer. JISRP; 9: 22–37.
6. Clayton MP. (2001) The effect of exhaled particles on RPE fit testing using the TSI Portacount plus respirator fit tester. PE/01/01. Bootle: HSE. Unpublished report.

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## What NIOSH Did

The purpose of the study was to find out if any of the QNFT methods could predict actual human exposure to a gaseous hazard. **Test subjects were fit tested and also exposed to Freon gas.** The respirators were equipped with filter cartridges that prevent Freon from passing through, therefore **any Freon detected in the test subject's bloodstream had to be due to inhalation and respirator leakage.**

Comparison of the different fit test methods was done by calculating a correlation coefficient based on the **measured fit factors and measured Freon blood concentrations. A coefficient of 1.0 means that there is a perfect match. Coefficients above 0.8 are usually considered to be very good and those below 0.5 are considered poor.** It is also necessary to compare coefficients to each other rather than simply looking at the value by itself because experimental uncertainties can shift the values, yet leave their relationship to each other intact.

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The NIOSH results for half mask respirators (Part II of NIOSH study)

The generated aerosol method (CHD) came out on top with a coefficient of 0.81. This confirmed it as the Gold Standard.

The PortaCount fit tester method (AA1) was a close second with a nearly identical coefficient of 0.78.

The CNP method yielded an unimpressive 0.36.

**NIOSH Results for Half Mask Respirators (Part II)**

QNFT Method	Correlation Coefficient ( $R^2$ )
Generated Aerosol (CHD)	0.81
PortaCount Fit Tester (AA1)	0.79
Controlled Negative Pressure (CNP)	0.36

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## What about Full-Face Respirators?

The poor results for the CNP method with half face respirators was corroborated in 2002 when NIOSH performed similar experiments using full-face respirators<sup>3</sup>. Those results cannot be directly compared to the half mask results because the coefficients were shifted due to some unknown reason. However, the relative differences between the coefficients confirm the sub-par performance of the CNP method discovered in Part II for half mask respirators. Relatively speaking, coefficients for the aerosol methods far exceeded the CNP coefficient. For full-face masks, the CHD and AA1 methods had coefficients of 0.09 and 0.11 respectively. The CNP coefficient was essentially zero ( $< 0.01$ ).

NIOSH Results for Full Face Mask Respirators

QNFT Method	Correlation Coefficient ( $R^2$ )
Generated Aerosol (CHD)	0.09
PortaCount Fit Tester (AA1)	0.11
Controlled Negative Pressure (CNP)	0